

Diameter and Base curve choice for the first contact lens

1. Measurement of the corneal diameter (HVID + 0.6 mm)* and K-readings
2. Determine the contact lens diameter \varnothing_T (use table below)
3. Determine of the Base curve $r_0 = r_{cfI} + BCf$ (use table below, r_{cfI} = flattest central K)

| | | Corneal diameter | | | | | | | |
|-----------------------|-------|-------------------------|-------|--------|-------|-------|-------|-------|-------|
| | | small | | medium | | | large | | |
| | | 11.10 | 11.30 | 11.50 | 11.70 | 11.90 | 12.10 | 12.30 | 12.50 |
| Contact lens diameter | 13.20 | 0.60 | 0.50 | 0.40 | | | | | |
| | 13.40 | 0.70 | 0.60 | 0.50 | 0.40 | | | | |
| | 13.60 | 0.80 | 0.70 | 0.60 | 0.50 | 0.40 | | | |
| | 13.80 | 0.90 | 0.80 | 0.70 | 0.60 | 0.50 | 0.40 | | |
| | 14.00 | 1.00 | 0.90 | 0.80 | 0.70 | 0.60 | 0.50 | 0.40 | |
| | 14.20 | 1.10 | 1.00 | 0.90 | 0.80 | 0.70 | 0.60 | 0.50 | 0.40 |
| | 14.40 | 1.20 | 1.10 | 1.00 | 0.90 | 0.80 | 0.70 | 0.60 | 0.50 |
| | 14.60 | | 1.20 | 1.10 | 1.00 | 0.90 | 0.80 | 0.70 | 0.60 |
| | 14.80 | | | 1.20 | 1.10 | 1.00 | 0.90 | 0.80 | 0.70 |
| | | Base Curve factor (BCf) | | | | | | | |

■ Orbis ($\varnothing_{\text{Cornea}} + 2.10 \text{ mm} / BCf = 0.60 \text{ mm}$)

■ Toris Bal – Torelis Bal – Borelis ($\varnothing_{\text{Cornea}} + 2.30 \text{ mm} / BCf = 0.70 \text{ mm}$)

■ Toris Int/Ext – Torelis Int/Ext ($\varnothing_{\text{Cornea}} + 2.50 \text{ mm} / BCf = 0.80 \text{ mm}$)

Example: Parameter for Toris Ballast:

Cornea parameters: $\varnothing_{\text{Cornea}} = 11.70 \text{ mm} / K_{\text{reading}} = 7.80 / 7.70 \text{ mm}$

- $\varnothing_T = 11.70 \text{ mm} + 2.30 \text{ mm} = 14.00 \text{ mm}$
- $r_0 = 7.80 \text{ mm} + 0.70 \text{ mm} = 8.50 \text{ mm}$

for 0.40 mm delta K, reduce 0.10 mm on r_0

Definitive 74: 0.10 mm steeper

* Information: 80% of the corneal curves are statistically between 11,3 and 12,1 mm.

Progress of the adaptation

1. Insert trial lens for a duration of between 30 minutes and 2 hours. Over refraction (you can use the autorefractometer for getting an idea of cyl/axis).
2. Biomicroscopy ($\times 10$ to 15) white light: observe the lens with patient looking straight ahead and during eye movement.
3. Mobility by eyelid movement, (Push up test).
4. Sag of the lens should be from 1 to 2 mm downwards.
5. Appearance of the front optic zone: tear film, hydration, lubrication, deposit.
6. Keratometry on the contact lens: (deformation of the mires).
7. Check for corneal and conjunctival staining with fluorescein after lens removal.
8. Order the definitive lens on the basis of the SN.