

HydroColor Pro an individually customized prosthetic lens with an exclusive white opaque back layer, available only from SwissLens

- a large selection of natural-looking standard colors or individually customized
- thin and comfortable contact lenses made of high-quality material (CTF 67)
- individual contact lens parameters
- iris and pupil diameters in different sizes
- hand dyed, with clear or black pupil
- consistently reproducible and colorfast

The colors do not affect the DK value of the contact lens.



HydroColor Pro - Selection

Our specially developed colour fan contains 30 different colors, enabling you to choose the optimal colour together with your customer on-site.

Diameter	iris	10.8 / 11.2 / 11.5 / 12.0 / 12.5 / 12.8
	clear pupil	2.7 / 3.3 / 3.8 / 4.2 / 4.8
	black pupil	2.6 / 3.2 / 3.7 / 4.2 / 4.8 / 5.2
	colour variants	30
Parameters	power	sph +/- 8.00 D, cyl -2.00 D
	Ø	13.00 to 16.00 mm
	ro	8.00 to 9.50 mm

Other parameters via Selection plus

HydroColor Pro - Creative

- individually selectable contact lens coloring
- order by color sample

To prevent color damage, especially fading of the opaque layer, please use only recommended solutions. The Acucare one step-T solution has been tested and is optimally compatible. Feel free to inquire about other recommended care products.

HydroColor Tint

We offer an additional option to order our soft contact lenses with a handling tint customized to your preferences. The following contact lens types are available:

Orbis, Toris, Borelis, Torelis, HydroMed, Relax

Transparent colors do not affect the DK value. Additionally, handling tints can temporarily change the eye color as a side effect.

Power

spherical:	+/- 20.00 D
toric:	sph +/- 20.00 D cyl. -4.00 D
Add:	+0.50 to +4.00 D
Material:	Contaflex CTF 67

We can accommodate additional corrections, special designs, and colors upon request.


ø Tint	11.0 / 11.3 / 11.8 12.2 / 12.6 / 13.0	2.6 / 3.2 / 3.7 / 4.2 / 4.8 / 5.2 6.2 / 7.2 / 8.0 / 9.0 / 10.0 / 11.0 11.3 / 11.8 / 12.2 / 12.6 / 13.0
ø clear pupil	2.6 / 3.2 / 3.6 4.2 / 4.7 / 5.8	-
Colours	1 2 3	1 2 3
blue		
aqua		
green		
brown		
red		
yellow		

HydroColor Intense

The HydroColor Intense is a light-reducing contact lens that can also be used for eye conditions such as achromatopsia.

Intense	Intense 1	Intense 2	Intense 3
brown			
grey			
brownS			

Caution: HydroColor Intense does not protect the conjunctiva.



Images are for illustration purposes only, final shade may vary.

HydroColor Filter

The HydroColor Filter is an edge-filter contact lens, serving as an alternative to glasses. It is available in various geometries, including spherical, toric, and multifocal.

- **edge filter contact lenses**

SLF 450	SLF 511	SLF 527	SLF 550

Filter selection should be based on Multilens eyeglass filters

- **blue light-reducing/blocking contact lens filters available upon request only.**

Parameters like HydroColor Tint, material CTF67